

The chemical structure of compound 1 is a complex molecule. It features a bicyclic core, specifically a decalin derivative. Substituents include a thioether group (-S-CH₂-CH₂-N(CH₃)CH₂CH₃) at position 1, a carboxylic acid group (-COOH) at position 2, and a hydroxyl group (-OH) at position 3. There are also methyl groups at positions 4, 5, and 6, and a cyclopropane ring fused to the bicyclic system at positions 7 and 8.

CCN(C)CCSC(=O)OC1C(C)C(C)C(C)C(C)C1C2C(C)C(C)C(C)C2OCCCC[N+](C)(C)Cc1ccccc1.[Br-]

Reference: PA/PH/81-108/T 0119 ANB

XXX:1724

COc1ccc(cc1)CNC(CO)c2cc(O)c(NC=O)cc2 and enantiomer
OC(=O)/C=C/C(=O)O 2 H₂O

Solvent mixture. Acetonitrile R. solution A (16:84 V/V).

The following chromatogram is published for information.




Figure 1724-1. Chromatogram obtained with reference solution (a) in the test for related substances.

Test solution. Dissolve 20.0 mg of the substance to be examined in the solvent mixture and dilute to 100.0 ml with the solvent mixture. This solution must be injected within 4 h from preparation, or stored protected from light at 4 °C for not more than 24 h.

Reference solution (a). Dissolve 5 mg of formoterol fumarate for system suitability CRS in the solvent mixture and dilute to 25.0 ml with the solvent mixture.

Reference solution (b). Dilute 1.0 ml of the test solution to 25.0 ml with the solvent mixture. Dilute 1.0 ml to 20.0 ml with the solvent mixture.

Column:

- size: $l = 0.15$ m, $\phi = 4.6$ mm;
- stationary phase: spherical octylsilyl silica gel for chromatography R3 (5 μ m)²⁰ with a pore size of 8 nm.

Mobility phase:

- mobile phase A: acetonitrile R1;
- mobile phase B: dissolve 3.73 g of sodium dihydrogen phosphate monohydrate R and 0.35 g of phosphoric acid R in water R and dilute to 1000 ml with the same solvent; the pH is 3.1 ± 0.1 .

Time (min)	Mobile phase A (per cent V/V)	Mobile phase B (per cent V/V)
0 - 10	25	84
10 - 37	15 → 70	84 → 30
37 - 40	70 → 15	30 → 84
40 - 65	15	84

Flow rate: 1.0 ml/min.

Detection: spectrophotometer at 214 nm.

Injection: 20 μ l.

Retention time: formoterol = about 12 min.

System suitability: reference solution (a):

- resolution: minimum 1.5 between the peaks due to impurity G and to impurity A.

Limits:

- correction factor: for the calculation of content, multiply the peak area of impurity A by 1.75;
- impurity A: not more than 1.5 times the area of the principal peak in the chromatogram obtained with reference solution (b) (0.3 per cent);
- impurity B, C, D or F: for each impurity, not more than the area of the principal peak in the chromatogram obtained with reference solution (b) (0.2 per cent); use the chromatogram obtained with reference solution (a) to identify the corresponding peaks;
- any other impurity: not more than 0.5 times the area of the principal peak in the chromatogram obtained with reference solution (b) (0.1 per cent);
- total: not more than 2.5 times the area of the principal peak in the chromatogram obtained with reference solution (b) (0.5 per cent);
- disregard limit: 0.25 times the area of the principal peak in the chromatogram obtained with reference solution (b) (0.05 per cent).

Impurity I. Liquid chromatography (2.2.29).

Test solution. Dissolve 5.0 mg of the substance to be examined in water R and dilute to 50.0 ml with the same solvent.

Reference solution (a). Dissolve 5.0 mg of formoterol impurity I CRS in water R and dilute to 50.0 ml with the same solvent.

Reference solution (b). Dilute 1.0 ml of the test solution to 20.0 ml with water R. Dilute 1.0 ml of this solution to 25.0 ml with water R.

(19) Karlson SDC, is suitable.

The following chromatogram is published for information.




Figure 1724.2. - Chromatogram obtained with reference solution (a) in the test for impurity I

Column:

- size: $l = 0.15$ m, $\phi = 4.6$ mm.

- stationary phase: octadecylsilyl vinyl polymer for chromatography *R₁₈*.

Mobile phase: mix 12 volumes of acetonitrile *R₁* with 88 volumes of a 5.3 g/l solution of tripotassium phosphate trihydrate *R* previously adjusted to pH 12.0 ± 0.1 with a 280 g/l solution of potassium hydroxide *R* or phosphoric acid *R*.

Flow rate: 0.5 ml/min.

Detection: spectrophotometer at 225 nm.

Injection: 20 μ l.

System suitability: reference solution (a):

- peak-to-valley ratio: minimum 3.5, where H_p = height above the baseline of the peak due to impurity I and H_v = height above the baseline of the lowest point of the curve separating this peak from the peak due to formoterol.

Limits:

- **Impurity I:** not more than the area of the principal peak in the chromatogram obtained with reference solution (b) (0.2 per cent).
- **disregard limit:** 0.25 times the area of the principal peak in the chromatogram obtained with reference solution (b) (0.05 per cent).

Water (2.5.12): 4.0 per cent to 5.0 per cent, determined on 50 mg.

ASSAY

Dissolve 0.350 g in 50 ml of anhydrous acetic acid *R*. Titrate with 0.1 M perchloric acid, determining the end-point potentiometrically (2.2.20).

1 ml of 0.1 M perchloric acid is equivalent to 40.24 mg of $C_{22}H_{23}N_3O_8$.


STORAGE

Protected from light.

IMPURITIES

Qualified impurities: A, B, C, D, E, F, I.

Other detectable impurities: G, H.



A. $R_1 = R_2 = R_4 = H$, $R_3 = CH_3$: 1-(3-amino-4-hydroxyphenyl)-2-[[2-(4-methoxyphenyl)-1-methylethyl]amino]ethanol,

B. $R_1 = CHO$, $R_2 = R_3 = R_4 = H$: *N*[[2-hydroxy-5-[(1*R*,3*S*)-1-hydroxy-2-[[2-(4-methoxyphenyl)ethyl]amino]ethyl]phenyl]formamide,

C. $R_1 = CO-CH_3$, $R_2 = R_4 = H$, $R_3 = CH_3$: *N*[[2-hydroxy-5-[(1-hydroxy-2-[[2-(4-methoxyphenyl)-1-methylethyl]amino]ethyl]phenyl]acetamide,

D. $R_1 = CHO$, $R_2 = R_3 = CH_3$, $R_4 = H$: *N*[[2-hydroxy-5-[(1-hydroxy-2-methyl[2-(4-methoxyphenyl)-1-methylethyl]amino]ethyl]phenyl]formamide,

E. $R_1 = CHO$, $R_2 = H$, $R_3 = R_4 = CH_3$: *N*[[2-hydroxy-5-[(1-hydroxy-2-[[2-(4-methoxy-3-methylphenyl)-1-methylethyl]amino]ethyl]phenyl]formamide,

(18) Available ODF-50 or Astec Polymer C₁₈ are suitable.

CN[C@H](C)Cc1ccc(OC)cc1 and enantiomer

O=C1C=CC(=C(C=C1)N)C2=CC=C(C=C2)C(O)CN(C)C3=CC=C(C=C3)C4=CC=C(C=C4)OC and enantiomerCOc1ccc(cc1)C[C@H](N)C[C@@H](O)c2ccc(NC=O)cc2 and enantiomer

Reagents

Tripotassium phosphate trihydrate. $K_3PO_4 \cdot 3H_2O$.
XXXXXXX. [22763-03-7].
White or almost white crystalline powder, freely soluble in water.


White or almost white crystalline powder, freely soluble in water.

Reference: PA/PH/Exp. 007490-19 AND JR

A new draft monograph on Metolazone, previously published in *Pharmeuropa* 11:3, is proposed.

XXXX:175T

Metolazonum



M.365.8

Draw the baseline between the transmission maxima at 813 cm^{-1} and at 649 cm^{-1} . Calculate R using the expression:

$$1 - \frac{T_1 - T_2}{T_3 - T_2}$$

T_1 = percentage transmittance of the baseline of the transmission minimum at 801 cm^{-1} .

Related substances. Liquid chromatography (2.2.29).

Test solution (a). Dissolve 320.0 mg of the substance to be examined in *dimethylformamide R* and dilute to 100.0 ml with the same solvent.